5 | Páfor støttebelte:

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1 | Scope of delivery

1 x Support belt »HEAVY«, with OUTLAST®

Read all of these instructions carefully!

2 | Intended use

By external pressure (compression), the abdominal and back muscles should be activated and supported. Supports the spine and promotes an upright posture.

- Disc problems
- Sciatica
- Lumbago
- Pain from overstrain

For private use.

3 | Contraindications

With the following symptoms, wearing the bandage is only recommended after consulting your doctor:

- Skin problems/injuries to the part of the body where the support is worn, especially if there are signs of inflammation, swelling, redness and overheating
- Sensation/movement disorders e.g., associated with diabetes
- open wounds and internal injuries
- allergic reactions
- Toddlers and people who can not be noticed when the bandage is uncomfortable

4 | Instructions for the safe use of the product

To date, there are no known side effects that impact the entire body. The product must be put on/used properly. The manufacturer is exempt from product liability if the product is used improperly.

- If the bandage is uncomfortable, remove it immediately!
- Wear bandage only by the hour, do not wear permanently!
- If your symptoms worsen or you notice any unusual changes, please contact your doctor immediately.

Do not use:

- while sleeping

5 | Applying support belt

Fig. 1: Hold the support belt with the hook and loop closure (inside) in the right hand and guide it around the back. Pull the two straps side-ways away from the body. Fig. 2: First fold in the left side and then the right one. Fig. 3: Then close the hook and loop fastener at any point on the bandage.

6 | Care instructions:

Close the hook and loop fastener to prevent damage on other clothes. The cleaning is done by household washing, according to the care instructions.

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