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GB | Support belt »EXTRA STRONG«, with OUTLAST®

1 | Scope of delivery
1 x Support belt »EXTRA STRONG«, with OUTLAST®
Read all of these instructions carefully!

2 | Intended use
By external pressure (compression), the abdominal and back muscles should be activated and supported. Supports the spine and promotes an upright posture.

- Disc problems
- Sciatica
- Lumbago
- Pain from overstrain

For private use.
3 | Contraindications
With the following symptoms, wearing the bandage is only recommended after consulting your doctor:

- Skin problems/injuries to the part of the body where the support is worn, especially if there are signs of inflammation, swelling, redness and overheating
- Sensation/movement disorders e.g., associated with diabetes
- open wounds and internal injuries
- allergic reactions
- Toddlers and people who can not be noticed when the bandage is uncomfortable

4 | Instructions for the safe use of the product
To date, there are no known side effects that impact the entire body. The product must be put on/used properly. The manufacturer is exempt from product liability if the product is used improperly.

- If the bandage is uncomfortable, remove it immediately!
- Wear bandage only by the hour, do not wear permanently!
- If your symptoms worsen or you notice any unusual changes, please contact your doctor immediately.

Do not use:
• while sleeping

5 | Applying support belt
Individually adjustable support via 2 additional elastic bands and 2 pull cords.

Fig. 1: First loosen the hook and loop fasteners of the two extra pull straps from the support belt. Hold the support belt with the hook and loop closure inside in the right hand and guide it around the back. Pull the two straps sideways away from the body. Fig. 2: First fold in the left side and then the right one. Then close the hook and loop fastener at any point on the bandage. Fig. 3+4: To attach the two extra strings proceed the same way.

6 | Care instructions:
Close the hook and loop fastener to prevent damage on other clothes. The cleaning is done by household washing, according to the care instructions.

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