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1 | Introduction

Thank you for deciding on this medical device.

The pelvic floor muscles are those muscle strands running from the pubic bone to the coccyx and forming a kind of loop around the genitals: around the female urethra, vagina and anus or around the male anus and penis base.

Above all for women a targeted training of these muscles is important for several reasons:

- A trained pelvic floor increases the pleasurable sensation and the orgasm capability during sex – the better the muscles are flexed and relaxed, the more joyful becomes the sex
- For the prevention of uterine descent women can become active themselves. Pelvic floor training is the ideal preparation for a (natural) birth. Above all pregnancy and birth can be responsible for the straining of the pelvic floor muscles. A targeted training helps to strengthen the muscles again which become strongly stretched during birth. Consequently, this helps enjoying sex soon again after birth
- Targeted pelvic floor training helps to prevent the muscles getting tired particularly during the menopause
- Trained pelvic floor muscles of women and men decrease or prevent incontinence, i.e. involuntary urination – for men, incontinence is often a problem above all after prostate operations

It is never too late to start training your pelvic floor. If starting with 20 or 70, it is just important to carry out your exercises regularly. Make the exercises part of your daily life.

2 | Important safety hints

Retain the manual for further use!


- To avoid a pinching of the skin in the thigh area it is recommended to use the Pelvic Muscle Trainer over clothes
- In case of any recent surgery, please consult your doctor before trying the Pelvic Muscle trainer. A physician should decide how long after a surgical procedure the muscles should be rested. Ask your doctor

3 | Main features



The Pelvic Muscle Trainer serves for a targeted workout of these muscle groups. The goal is to challenge your muscles a little more each day, not to tire them out or create pain.

3.1 | Uses

1. Warm up before and wear comfortable clothes while using Pelvic Muscle Trainer.
2. Place Pelvic Muscle Trainer between thighs, few inches below pelvic floor. **DO NOT** place against pelvic floor.
3. Center the Pelvic Muscle Trainer with Wing flaps touching front or back of thighs depending which is more comfortable
4. Maintain proper posture. Stand up straight with feet close together, knees slightly bent.
5. While keeping knees, hips and shoulders parallel, close Wing Enclosure Area by squeezing Wing Flaps with upper thighs, buttocks and pelvic floor muscles.
6. For beginners, squeeze and hold the Pelvic Muscle Trainer for 3 seconds and then rest 7 seconds. Try to perform 6 repetitions. You may find it difficult to close the Wing Enclosure Area at first. This is OK because you are still activating the PC muscles even when you get the Wing Flaps halfway closed.
7. Gradually increase usage and squeeze hold time. Start by using the Pelvic Muscle Trainer three times a week and increase to twice daily. Give with enough time for the muscles to recuperate. **AVOID STRAINING THE MUSCLES.**
8. Experienced users should be able to perform 10 repetitions. Holding 25 seconds resting 5 seconds
9. As long as you maintain posture, exercise regularly and keep the Pelvic Muscle Trainer positioned properly the strength will come automatically.

Always make sure you close the Flex Flaps as far as you can with each repetition. The training is efficient and has no side effects. Effects are individually variable depending on the regularity of use.

Other applications for training the inner thigh muscles or the chest muscles:

4 | Cleaning and storing

- Store in a clean, dry place
- Do not place or store in direct sunlight
- After each use wipe the device with a dry or damp cloth

5 | Delivery Contents

- Pelvic Muscle Trainer
- Instruction manual

6 | Technical Data

- Model: SKE2
- Material: Polypropylen
- Weight: 205 gr
- Dimensions: 20 cm x 23 cm x 8 cm

This product may be safely disposed of in normal home trash

7 | Warranty

The guarantee lasts for 24 months from the date of purchase for manufacturers' error in material and quality. If it is used for another purpose or operated incorrectly, no liability for possible damages can be accepted. This guarantee is only valid if the instructions in the "instructions for use" have been followed.

8 | Returns

Dear Customer,

In case your product does not function as you expected please do not return immediately. Perhaps our Hotline (+49-69-954061-24) can be of assistance. If you are still not satisfied, you have the possibility to return products free of charge within the BRD. We require your address. This you may provide to the above listed telephonic hotline or by email to: technik@hydas.de. Prepare the package for shipping and include the receipt and description of the problem. Please be aware of the fact that collect packages may not be accepted.

Your Hydas customer service team

Address for the guarantee and service:

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Hydas GmbH
c/o atrikom fulfillment GmbH • Carl-Zeiss-Straße 10
Hydas GmbH & Co.KG
Hydas GmbH & Co. KG
c/o atrikom fulfillment GmbH